

# Traditional Wound Management



## Cosmopor®

Adhesive wound dressing



- Soft and gentle to the skin – can be used for patients with sensitive skin
- Conforms well to body contours and remains in place
- Rounded corners help reduce the need to change the dressing
- Good absorbency and cushioning effect
- Fluid repellent backing

## Peha-Haft®

Conforming cohesive bandage



- Conforming elastic bandage that sticks to itself and not the skin, hair or clothes
- Very kind to the skin
- Does not contain natural rubber latex, chlorine or acrylate
- No need for additional fixation
- High absorption and air permeable

## Omnifix®

Nonwoven retention/abrasion tape



- Skin-friendly adhesive
- Especially soft and conformable
- Highly permeable to air and water vapor
- Helps reduce potential for maceration
- Provides reliable adherence
- Self-adhesive, nonwoven

## Sterilux® Bulky

Premium bulky bandage



- Made of 100% highly absorbent gauze
- Air permeable
- Made without optical brightener or chlorine
- Sterile
- Does not contain natural rubber latex


## Flexicon®

Premium conforming stretch bandage



- Soft, polyester fabric for added comfort
- Holds dressings securely in place without slipping
- Unique structure allows bandage to cling to itself
- Conforms to difficult-to-dress body contours
- Finished edges minimize lint and unraveling
- Blue application guide helps ensure consistent wrapping

## Shapes:



## Indications:

- Protection of post-surgical wounds
- Wounds with low to moderate level of exudate
- Protection of wounds following wound cleaning

## Sizes:

- 1" x 4.5 yards
- 1.5" x 4.5 yards
- 2.25" x 4.5 yards
- 3" x 4.5 yards
- 4" x 4.5 yards



## Indications:

- Safe dressing retention for use on joints and conical or rounded parts of the body
- Fixation of cannulae, padding material, etc.

## Sizes:

- 2" x 2 yards
- 4" x 2 yards
- 2" x 10 yards
- 4" x 10 yards
- 6" x 10 yards



## Indications:

- Dressing retention
- Ideal for use on joints, round and angular parts of the body
- Securement of catheters, measuring instruments, probes and similar devices.

## Sizes:

- Sterile  
4.5" x 4.1 yards
- Non-Sterile  
4.5" x 4.1 yards




## Indications:

- Primary treatment of dirty, infected and highly exudative wounds
- Secondary dressing
- Universal wound dressing

## Sizes:

- | Sterile        | Non-Sterile    |
|----------------|----------------|
| 1" x 4.1 yards | 1" x 4.1"      |
| 2" x 4.1 yards | 2" x 4.1 yards |
| 3" x 4.1 yards | 3" x 4.1 yards |
| 4" x 4.1 yards | 4" x 4.1 yards |
| 6" x 4.1 yards | 6" x 4.1 yards |



## Indications:

- Applications that require controlled compression
- Securing dressings
- Treating burns





## Conco®

Conforming stretch bandage

- Secures dressings in place with minimal risk of slipping or constriction
- Clings to itself when overlapped
- Conforms to body contours
- Finished edges minimize lint and unraveling
- Does not contain natural rubber latex

## Sizes:

| Sterile        | Non-Sterile    |
|----------------|----------------|
| 1" x 4.1 yards | 2" x 4.1 yards |
| 2" x 4.1 yards | 3" x 4.1 yards |
| 3" x 4.1 yards | 4" x 4.1 yards |
| 4" x 4.1 yards | 6" x 4.1 yards |
| 6" x 4.1 yards |                |



## Indications:

- Applications that require controlled compression
- Effective for securing dressings
- For use in treating burns



## Sorbalex® ABD

Absorbent abdominal pad

- Blue guide line to ensure proper application
- Absorbent cellulose inner layer core
- Soft, nonwoven outer layer for patient comfort
- Does not contain natural rubber latex

## Shapes:

|           |
|-----------|
| 5" x 9"   |
| 7.5" x 8" |
| 8" x 10"  |



## Indications:

- Primary or secondary dressing for heavily draining acute or chronic wounds
- Secure with Ominfix® latex-free retention tape



## Sorbalex® Non-Adherent

Absorbent non-adhering dressing

- Bonded on both sides by a smooth, perforated polyester film outer surface
- Minimizes risk of disrupting tissue healing by not sticking to the wound
- Allows fluids to pass through dressing
- Cuts to fit various wound sizes with virtually no loose fibers
- Absorbent inner layer of polyester/ rayon material

## Shapes:

|         |
|---------|
| 2" x 3" |
| 3" x 4" |
| 3" x 8" |



## Indications:

- Primary dressing for wounds such as burns and early stage skin ulcers
- Primary dressing for lightly draining acute or chronic wounds



## Sterilux®

Premium gauze sponges

- 100% highly absorbent gauze
- Can be unfolded several times to obtain a larger size without fraying
- Made with pure, clean cotton without optical brightener or chlorine
- Offered in different product sizes, layers and threads
- Sterile

## Shapes:

Sterile      Non-Sterile

|         |         |
|---------|---------|
| 2" x 2" | 2" x 2" |
| 3" x 3" | 3" x 3" |
| 4" x 4" | 4" x 4" |

## Indications:

- Universal wound dressing
- Wound cleansing or as a wound cover
- Dirty, infected and high level of exudate
- Secondary dressing

## Thickness:

8 ply

12 ply






**Econolux®**  
Gauze sponges

- 100% absorbent gauze
- Cut edges are on the inside of the sponge
- Offered in different product sizes, layer and threads
- Sterile

**Shapes:**

Sterile                  Non-Sterile



**Indications:**

- General wound treatment
- Wound cleansing or as a wound cover
- Dirty, infected and high level of exudate
- Secondary dressing




**Sorbalux® All-Purpose**

All-purpose drain and I.V. cover sponges

- Soft polyester/ rayon nonwoven materials
- Absorbent and fast wicking
- Pre-cut slit drain and I.V. sponges
- 4 - and 6-ply options
- Easy to apply


**Shapes:**

Sterile




**Shapes:**


Non-Sterile



Drain sponge  
4" x 4"



I.V. sponge  
2" x 2"



**Indications:**

- Wound cleansing and prepping, applying ointment wiping needles, cleaning slides, etc.
- Secure with Ominfix® latex-free retention tape
- All-Purpose Sponge is ideal as a primary or secondary wound dressing
- Drain Sponge is ideal as sterile protection around chest, tracheotomy, feeding and gastric tubes
- I.V. Non-Woven Sponge is ideal as sterile protection around I.V., mid-line and central line catheters



**Flex-Band®**


Fabric adhesive bandage

- Skin friendly
- Special shapes to allow for optimal movement
- Does not adhere to the wound bed
- Flexible and snug fitting
- Good adhesion
- Sterile
- Does not contain natural rubber latex

**Indications:**

- Protection of minor injuries
- Protect wounds including the knuckle or finger tip

**Shapes:**



# Ordering information

| Order no. | Product name                 | Size            | Pcs/<br>Pack | Pcks/<br>Case |
|-----------|------------------------------|-----------------|--------------|---------------|
| 22100000  | Flexicon® Non-Sterile        | 1" x 4.1 yds.   | 24           | 4             |
| 22200000  | Flexicon® Non-Sterile        | 2" x 4.1 yds.   | 12           | 8             |
| 22300000  | Flexicon® Non-Sterile        | 3" x 4.1 yds.   | 12           | 8             |
| 22400000  | Flexicon® Non-Sterile        | 4" x 4.1 yds.   | 12           | 8             |
| 22600000  | Flexicon® Non-Sterile        | 6" x 4.1 yds.   | 6            | 8             |
| 19100000  | Flexicon® Sterile            | 1" x 4.1 yds.   | 12           | 8             |
| 19200000  | Flexicon® Sterile            | 2" x 4.1 yds.   | 12           | 8             |
| 19300000  | Flexicon® Sterile            | 3" x 4.1 yds.   | 12           | 8             |
| 19400000  | Flexicon® Sterile            | 4" x 4.1 yds.   | 12           | 8             |
| 19600000  | Flexicon® Sterile            | 6" x 4.1 yds.   | 12           | 4             |
| 18200000  | Flexicon® Clean Wrap         | 2" x 4.1 yds.   | 20           | 5             |
| 18300000  | Flexicon® Clean Wrap         | 3" x 4.1 yds.   | 20           | 5             |
| 18400000  | Flexicon® Clean Wrap         | 4" x 4.1 yds.   | 20           | 5             |
| 18600000  | Flexicon® Clean Wrap         | 6" x 4.1 yds.   | 20           | 3             |
| 80100000  | Conco® Non-Sterile           | 1" x 4.1 yds.   | 24           | 4             |
| 80200000  | Conco® Non-Sterile           | 2" x 4.1 yds.   | 12           | 8             |
| 80300000  | Conco® Non-Sterile           | 3" x 4.1 yds.   | 12           | 8             |
| 80400000  | Conco® Non-Sterile           | 4" x 4.1 yds.   | 12           | 8             |
| 80600000  | Conco® Non-Sterile           | 6" x 4.1 yds.   | 6            | 8             |
| 81200000  | Conco® Sterile               | 2" x 4.1 yds.   | 12           | 8             |
| 81300000  | Conco® Sterile               | 3" x 4.1 yds.   | 12           | 8             |
| 81400000  | Conco® Sterile               | 4" x 4.1 yds.   | 12           | 8             |
| 81600000  | Conco® Sterile               | 6" x 4.1 yds.   | 12           | 4             |
| 83400000  | Sterilux® Bulky Non-Sterile  | 4.5" x 4.1 yds. | 1            | 100           |
| 83500000  | Sterilux® Bulky Sterile      | 4.5" x 4.1 yds. | 1            | 100           |
| 56810000  | Sterilux® Non-Sterile 8 ply  | 2" x 2"         | 200          | 25            |
| 56820000  | Sterilux® Non-Sterile 12 ply | 2" x 2"         | 200          | 40            |
| 56840000  | Sterilux® Non-Sterile 12 ply | 3" x 3"         | 200          | 20            |
| 56850000  | Sterilux® Non-Sterile 8 ply  | 4" x 4"         | 200          | 20            |
| 56860000  | Sterilux® Non-Sterile 12 ply | 4" x 4"         | 200          | 10            |
| 56880000  | Sterilux® Sterile 8 ply      | 2" x 2"         | 100          | 30            |
| 56890000  | Sterilux® Sterile 12 ply     | 3" x 3"         | 80           | 30            |
| 56900000  | Sterilux® Sterile 8 ply      | 4" x 4"         | 100          | 12            |
| 56910000  | Sterilux® Sterile 12 ply     | 4" x 4"         | 50           | 24            |
| 416101    | Sterilux® Sterile 8 ply      | 4" x 4"         | 10           | 128           |
| 416814    | Econolux® Non-Sterile        | 2" x 2"         | 200          | 25            |
| 416815    | Econolux® Non-Sterile        | 2" x 2"         | 200          | 40            |
| 416816    | Econolux® Non-Sterile 12 ply | 4" x 4"         | 200          | 20            |
| 416817    | Econolux® Non-Sterile 16 ply | 4" x 4"         | 200          | 10            |
| 416818    | Econolux®                    | 4" x 4"         | 200          | 10            |
| 416103    | Econolux® Sterile (2s)       | 2" x 2"         | 100          | 30            |
| 416104    | Econolux® Sterile (2s)       | 4" x 4"         | 100          | 12            |
| 416105    | Econolux® Sterile (2s)       | 4" x 4"         | 50           | 24            |

| Order no. | Product name               | Size                     | Pcs/<br>Pack | Pcks/<br>Case |
|-----------|----------------------------|--------------------------|--------------|---------------|
| 48700000  | Sorbalux® ABD              | 5" x 9"                  | 25           | 16            |
| 48710000  | Sorbalux® ABD              | 7.5" x 8"                | 20           | 12            |
| 48720000  | Sorbalux® ABD              | 8" x 10"                 | 16           | 20            |
| 48800000  | Sorbalux®                  | 4" x 4" Drain            | 100          | 6             |
| 48810000  | Sorbalux®                  | 2" x 2" I.V.             | 140          | 10            |
| 48830000  | Sorbalux® AP               | 2" x 2"                  | 200          | 40            |
| 48840000  | Sorbalux® AP               | 3" x 3"                  | 200          | 20            |
| 48850000  | Sorbalux® AP               | 4" x 4"                  | 200          | 10            |
| 48860000  | Sorbalux® Sterile (2s)     | 2" x 2"                  | 100          | 30            |
| 48870000  | Sorbalux® Sterile (2s)     | 3" x 3"                  | 100          | 24            |
| 48880000  | Sorbalux® Sterile (2s)     | 4" x 4"                  | 100          | 12            |
| 48890000  | Sorbalux® NA* Sterile (2s) | 2" x 3"                  | 100          | 24            |
| 48900000  | Sorbalux® NA* Sterile (2s) | 3" x 4"                  | 100          | 24            |
| 48910000  | Sorbalux® NA* Sterile (2s) | 3" x 8"                  | 50           | 12            |
| 932452    | Peha-haft®                 | 1" x 4.5 yds.            |              | 8             |
| 932441    | Peha-haft®                 | 1.5" x 4.5 yds.          |              | 90            |
| 932442    | Peha-haft®                 | 2.25" x 4.5 yds.         |              | 90            |
| 932443    | Peha-haft®                 | 3" x 4.5 yds.            |              | 120           |
| 932444    | Peha-haft®                 | 4" x 4.5 yds.            |              | 135           |
| 46100000  | Flex-Band®                 | 3" x 3" 4-Wing           | 50           | 24            |
| 46110000  | Flex-Band®                 | 7/8" Round Spot          | 100          | 24            |
| 46130000  | Flex-Band®                 | 2" x 3" Patch            | 50           | 24            |
| 46140000  | Flex-Band®                 | 1 1/2" x 3" Knuckle      | 100          | 24            |
| 46150000  | Flex-Band®                 | 1 1/4" x 1" Oval         | 100          | 24            |
| 46160000  | Flex-Band®                 | 3/4" x 3" Strip S        | 100          | 24            |
| 46180000  | Flex-Band®                 | 1" x 3" Strip L          | 100          | 24            |
| 46200000  | Flex-Band®                 | 2" x 1 3/4" Fingertip S  | 100          | 24            |
| 46210000  | Flex-Band®                 | 2" x 2 1/2" Fingertip L  | 50           | 24            |
| 36200000  | Omnifix®                   | 2" x 2 yds. unstretched  |              | 24            |
| 36400000  | Omnifix®                   | 4" x 2 yds. unstretched  |              | 24            |
| 36210000  | Omnifix®                   | 2" x 10 yds. unstretched |              | 24            |
| 36410000  | Omnifix®                   | 4" x 10 yds. unstretched |              | 12            |
| 36610000  | Omnifix®                   | 6" x 10 yds. unstretched |              | 12            |
| 900800    | Cosmopor® Sterile          | 2" x 2.8"                | 50           | 20            |
| 900806    | Cosmopor® Sterile          | 3.2" x 4"                | 25           | 22            |
| 900808    | Cosmopor® Sterile          | 3.2" x 6"                | 25           | 14            |
| 900812    | Cosmopor® Sterile          | 4" x 8"                  | 25           | 12            |
| 900814    | Cosmopor® Sterile          | 4" x 10"                 | 25           | 8             |
| 900816    | Cosmopor® Sterile          | 4" x 14"                 | 25           | 8             |
| 900820    | Cosmopor® Sterile          | 4" x 4"                  | 25           | 8             |
| 900823    | Cosmopor® Sterile          | 6" x 6"                  | 25           | 8             |

NA = Non-Adherent  
AP = All Purpose

To learn more about our products or to place an order contact  
HARTMANN at 1-800-243-2294 or visit [hartmannusa.com](http://hartmannusa.com).

HARTMANN USA, Inc.  
Rock Hill, SC 29730  
1-800-243-2294

[hartmannusa.com](http://hartmannusa.com)

XLIT 2699 Rev. 1 (01/17)



Going further  
for health