






RespoSorb® Super



Spec.-No.:	D 6.1203
Department:	DEV-PHF
Date:	2018-04-05

Schematic view:



Sterilization

Ethylen Oxide (EO) according to NF EN ISO 11135-1

5. Product Requirements

Absorbency Minimum 140g/100cm<sup>2</sup> of absorbent core surface, EN 13726-1 (chap 3.2)

Dimensions

Format (cmxcm)	10x10	10x20	15x20	20x25	20x40
Absorbent pad length (cm)	10	20	20	25	40
Absorbent pad width (cm)	10	10	15	20	20
Tolerance (mm)	5	5	5	5	5

6. Labelling

Lot-No. With 9+6 digit code:

x	XXX	XX	XX	X		XX	X	X	XX
Year	Internal order number	Charge of sterilization	plant	Control digit		week	day	shift	machine

Expiry date: Year Month

Example: 2018 06

Shelf life: 3 Years



# Technical Data Sheet

## RespoSorb® Super



Spec.-No.:	D 6.1203
Department:	DEV-PHF
Date:	2018-04-05

### 7. Packaging

RespoSorb® Super is individually packed in paper-paper peel-pouches (all sizes excepted 20x40) or in paper-film blister (size 20x40). These packaging are placed in folding boxes and transport carton according to DIN, sealed with adhesive tapes and packed onto Europallet.

Date: 05/04/2018

PAUL HARTMANN AG  
Product Development – PHF

Hélène Jouët-Pastré